




Ophthalmic Surgery



Ophthalmic Surgery

Dear Customer's

Thank you for choosing **Sadqain Enterprises** and taking the time to look at our newest Ophthalmic Instruments catalogue.

Since we have been working successfully in the medical field. We try to serve you better and make things easier for you. Please make best use of this catalogue for ordering products and services required.

We constantly try to upgrade the quality of instruments and equipments. We need your unending suggestions intended for the improvements of quality and services of our products. Your suggestions are necessary, to serve you better - pleasantly - efficiently and reliably.

Your suggestions to improve the contents of this catalogue are most welcome. C.E.O.
Sadqain Enterprises

Forceps

2 - 7

Scissors

8 - 9

Speculums

10 - 11

Choppers & Rotators

12 - 13

Markers

14 - 15

Vitreo Retinal Instruments

17

Curettes

18 - 19

Dialtors

19

Calipers

20

Hooks & Fixation Rings

21

Needle Holders

23


T= Available in Titanium

Moorields Forceps



SE-112

Blunt tipped forceps for manipulating the conjunctiva.




Small Handle Straight Tying Forceps



SE-119

Straight shafts with 4.5mm tying platforms, 85mm long.




Small Handle Colibri Tying Forceps



SE-120

1 into 2 teeth, 0.12mm with tying platforms, 85mm long.




Straight Toothed Forceps



SE-121

1 into 2 teeth, with 5.0mm tying platforms, 100mm long.




Colibri Style Toothed Forceps




SE-122

1 into 2 teeth, 0.12mm with tying platforms, 105mm long.



T= Available in Titanium


Short Handle Capsularhexis Forceps



SE-130

45° angled straight shafts. Tip to angled length 14.0mm. 85mm long.


Capsularhexis Forceps Straight



SE-131

45° angled straight shafts. Tip to angled length 14.0mm.


Capsularhexis Forceps Curved



MS132

45° angled gently curved shafts. Tip to angled length 14.0mm, 105mm long.


X Action Capsularhexis Forceps



SE-133

45° angled gently curved shafts. Tip to angled length 12.0mm.

Round Bodied Capsularhexis Forceps




SE-134


45° angled straight shafts. Tip to angled length 14.0mm.

T= Available in Titanium

Utrata Style Non-Toothed Capsularhexis Forceps




Swan Neck Type Plain Forceps



Straight Notched Forceps




0.25mm notched, with tying platforms, 115mm long.



Curved Notched Forceps




0.25mm notched, with tying platforms, 115mm long.



Lens Folder




Folding 5-6mm sized optics of intra-ocular lenses (IOL).
Paddle style jaws with handle stop to prevent over-folding.



T= Available in Titanium


Lens Implantation Forceps



SE-152

To implant the folded intra-ocular implant into the capsular bag.



Locking Lens Introducer



SE-153

Standard cross action with lock.


St Martins Forceps 1in 2 teeth



SE-154

1 into 2 teeth, with 5.0mm tying platforms, 85mm long.


Hoskins Forceps



SE-155

0.25mm notched, with tying platforms, 85mm long.

Bishops Forceps




SE-156

T= Available in Titanium

Epilation-Cilia Forceps



SE-158




Bishop-Harmon Tissue Forceps



SE-165

1 into 2 teeth, Delicate, 85mm long.




Straight Tying Forceps



SE-171

Straight shafts with 4.5mm tying platforms, 105mm long.




McPherson Tying Forceps



SE-172

Angled 45° shafts with 4.5mm tying platforms for tying knots and also manipulation of intra-ocular lens, 100mm long.




Kelmans Forceps




SE-173

Angled shafts with tying platforms.




T= Available in Titanium



Girard straight tying forceps

SE-174


Round handle, 105mm long.



Jewelers forceps

SE-175


0.1x0.06mm tips, Straight, 110mm long.



Scleral plug forceps


SE-177

118mm long.



Mosquito Type Forceps 12.5cm

12.5cm general holding forceps.




Curved Mosquito Forceps

T= Available in Titanium


Westcott Scissors




Ong Type Scissors




Vannas Scissors




Castroviejo Corneal Section Scissors (Right)



Castroviejo Corneal Section Scissors (Left)






General Scissors Sharp Pointed



Sharp pointed 9.5cm scissors.



General Scissors Blunt End




Blunt Ended 9cm scissors.



General Scissors Long Blunt End



Blunt ended 13cm scissors.




Steven Scissors Curved




T= Available in Titanium

Closed loop Kratz Barraquer wire speculum.




Closed loop Kratz Barraquer wire speculum.

SE-109



Kratz Type Paediatric 10mm closed wire Speculum



10mm blades.

SE-110




Kratz Barraquer Wire Speculum




Speculum of steel wire with 13.5mm open blades, 52mm long.

SE-111




Temporal Speculum




Speculum of steel wire, 15mm blades.

SE-112




Closed Loop Kratz Speculum



Closed loop speculum of steel wire.

SE-113



T= Available in Titanium

Sprung Speculum Paediatric 6mm blade



SE-114



6mm blades.

Sprung Speculum Neonate 4mm blade



SE-115




4mm blades.

Pierce Type Adjustable Speculum



SE-116



Slide adjustable speculum.

Bladed Adjustable Speculum




SE-110



Screw adjustable speculum with solid blades.

Adjustable Speculum



SE-111



Screw adjustable speculum, wire type.

T= Available in Titanium

Haeliger phaco cleaver



Y Nucleus Rotator



Claw Phaco Chopper




Phaco Chopper



Koch Type Chopper




T= Available in Titanium



Green Chopper

SE-145

Angled nucleus chopper.



Sinskey Lens Dialer / Sinskey Hook

SE-146


For manipulating the intra-ocular lens.



Drysdale Nucleus Spatula

SE-147

Paddle shaped tip.



Mushroom Nucleus Manipulator

SE-148


For iris retraction.

T= Available in Titanium

Lasik Corneal Marker




SE-121




Two blunt ins angled downwards and separated at 45° angle - to mark the cornea.

Flap Elevator



SE-131




Blunt tipped 1mm spatula for manipulating the lasik corneal lap.

Universal Flap Lifter




SE-132




Bevelled tip suitable for manipulating femto lasik laps.

Intra-Lasik Flap Separator



SE-134



Round, sharp tipped separator.

T= Available in Titanium


CK Lasik Marker (target center)



SE-140

Target center to facilitate accurate marking of the cornea.


Lasek Alcohol Well 9mm diameter



SE-141

Liquid holding well for alcohol application during surgery.


Lasek Alcohol Well 8mm diameter



SE-142

Liquid holding well for alcohol application during surgery.

Lasek Alcohol Well 7mm diameter



SE-143

Liquid holding well for alcohol application during surgery.


T= Available in Titanium

Epi-peeler



SE-151

For manipulating the epithelium in lasek procedures.



Hoe Type Epi-peeler (D Shaped)



SE-160

For manipulating the epithelium in lasek procedures.




Hoe Type Epi-peeler (Square Shaped)



SE-161

For manipulating the epithelium in lasek procedures.




Lasik Flap Forcep



SE-171

For manipulating the flap in lasik procedures.




T= Available in Titanium

Vitreo Retinal Instruments



Squeeze handle instruments for Vitreo Retinal surgery



Side Gripping Forceps
Smooth Round Jaws,
2.4mm.

SE-1501




Foreign Body Forceps
with cap jaws,
straight.

SE-1502




Gripping Forceps,
2.5mm Straight.

SE-1505




Gripping Forceps,
2.0mm Smooth
Round Tips.

SE-1506



Serrated, Straight.

SE-1507




Straight, and Gripping
Jaws.

SE-1509




Scissor, Curved,
(R12mm).

SE-1153




Scissor, Angled, 45°.

SE-1154



Scissor, Vertical, 40°.

SE-1155



Scissor, Vertical, 80°.

SE-1156

T= Available in Titanium

Chalazion Clamp Small



SE-101

11mm diameter, round plate.



Chalazion Clamp Medium



SE-102


Oval plate, 23x16mm.



Eye Lid Plate



SE-117




Curette Small



SE-101

Chalazion curette - 1mm diameter.



Curette Medium




SE-102

Chalazion curette - 1.5mm diameter.



T= Available in Titanium

Curette Large




SE-103



Chalazion curette - 2mm diameter.


Punctal Dilator Sharp



SE-101




Punctal Dilator Blunt



SE-102



Trabeculectomy Punch



SE-103



Bullet shaped tip for easy insertion into sclero-corneal tunnel.

Desmarres Retractor




SE-120




T= Available in Titanium

Braunstein Fixed Caliper



Braunstein Fixed Caliper - Fine




Castroviejo Adjustable Caliper




Scleral Marker



Jafe Type Lid Retractor Top Lid 14mm blades(x2)




T= Available in Titanium



SE-109


Jameson Strabismus Hook



SE-110

5mm hook.


Jameson Hook Extended



SE-111

8mm hook.


Fine Thornton Fixation Ring Small 13mm



SE-112

Swivel head 13mm/11mm cut out.

Fine Thornton Fixation Ring Standard 16mm



SE-113

Swivel head 16mm/14mm cut out.

T= Available in Titanium

Freer Elevator




Heveston muscle hook



Heveston muscle hook




T= Available in Titanium



SE-161

Needle Holder

Curved, 6mm delicate jaw.



SE-162

Castroviejo Needle Holder with lock

Ophthalmic Surgery



House No 344 Kotli Loharan East Sialkot

Cell: +92 322-4007797

Email: info@sadqainenterprises.pk

web: www.sadqainenterprises.pk

Ophthalmic Surgery